
# X-Ray Photoelectron Spectroscopy

## Samples handling

Surface Analysis | [pardis.simon@univ-lille.fr](mailto:pardis.simon@univ-lille.fr)

- **Identify your need of an XPS Analysis**  
Is XPS appropriate for what you are looking for?
- **Prepare an analysis form**  
Download it on the website
- **Come and discuss with us!**
- **Selection of samples of interest**  
The less the better!

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
Analyse de surface

Request date:  
 Organism:  
 Full name:  
 E-mail:  
 Phone number:  
 Supervisor name:  
 Billing address (please specify the contract name):

Supervisor Signature:

CONTEXT OF THE STUDY:

DETAILED NATURE OF SAMPLES \* (number, atomic composition, synthesis method, dimensions [length, width, thickness]):

ANALYSIS REQUESTED:

XPS – information required \*

- **Ultra-high vacuum compatible samples**
- Samples can be all kind of solids, thin films, powders, fibers, porous solids, foams...; but also ionic liquids with low vapor pressure (below  $10^{-10}$  mbar)...
- Conductors, insulators, ceramics, alloys, polymers, organic materials... can be analyzed by XPS

## Two rules:

- **Avoid contamination** as XPS is surface sensitive!

During each step: handling, sample storage, exposure to the atmosphere, preparation for analysis

- **Maximize electrical contact with analyzer** to avoid charging effect

Properly mounting the samples for analysis to create a constant surface potential

## Sample handling, preparation and mounting for XPS and other surface analytical techniques

Cite as: J. Vac. Sci. Technol. A 38, 063202 (2020); doi:10.1116/6.0000421

Submitted: 30 June 2020 - Accepted: 20 August 2020

Published Online: 11 September 2020

Fred A. Stevie,<sup>1,4\*</sup> Roberto Garcia,<sup>1</sup> Jeffrey Shallenberger,<sup>2</sup> John G. Newman,<sup>3</sup> and Carrie L. Donley<sup>4</sup>

### AFFILIATIONS

<sup>1</sup>Analytical Instrumentation Facility, North Carolina State University, Raleigh, North Carolina 27695

<sup>2</sup>Materials Characterization Laboratory, The Pennsylvania State University, University Park, Pennsylvania 16802

<sup>3</sup>Physical Electronics, Chanhassen, Minnesota 55317

<sup>4</sup>Chapel Hill Analytical and Nanofabrication Laboratory, University of North Carolina, Chapel Hill, North Carolina 27599

**Note:** This paper is part of the Special Topic Collection on Reproducibility Challenges and Solutions.

\*Electronic mail: fred.stevie@ncsu.edu


## Avoid contamination

- Store the samples in **small glass** vial or wrap them in aluminum foil – **avoid plastic**, high risk of carbon or siloxane contamination
- For sample sensitive to air or humidity, store in Ar filled air tight container
- ➔ Possibility to introduce the samples in the UHV chamber with a glove bag filled with N<sub>2</sub>.
- Samples **must be dry** – analysis is performed under UHV

## Avoid contamination

- Clean all the tools and sample holder with **acetone first then ethanol or isopropanol**
- Use **fiber-free** tissues
- **NEVER manipulate the samples with bare hands**, use nitrile or powder free/latex free gloves and don't touch sample surface with gloves or tweezers either, manipulate the edges
- For solids cut to the appropriate size with **clean, lubricants free tools** and protect the surface
- Blow the contamination particles **with dry air or N<sub>2</sub>**.

## Sample dimension

**Kratos Axis Ultra DLD Instrument** - Two types of sample holder adapted for solids and/or powders

- For solids cut to the appropriate size:

**Width, length < 4 cm ; Thickness < 4 mm**

- For powders or ionic liquids :

**A volume of ca 0.1 ml is needed**

# Maximize electrical contact

## Powders

## Solids

- Compressed in specific sample holder (compatible with ionic liquids)
- Double sided conductive tape (Copper or carbon)



- Metallic clips or metallic fixation plates with screw to fix the samples on the sample holder
- Double sided conductive tape (Copper or carbon)

Ionic liquids

Free powders

Taped powders or solids

Fixed solids

Adapted for measurements with ToF-SIMS/LEIS/High T° gas Cell

Adapted for angle resolved measurement